



JUVENT
Health

Indications For Use

The Juvent 1000N Platform Is Indicated For Maintaining and/or enhancing muscle strength, function and postural stability. Juvent's 25 year research history has shown that the technology:

- Restores dynamic action of the calf muscle fibers
- Strengthens & increases muscle mass
- Stimulates blood & lymphatic flow
- Improves balance, postural reflexes and coordination

Medical Device Precautions

There have been no reports of adverse effects with the use of the Juvent 1000N DMT Platform to date. This includes thousands of treatments in children, adolescents, and adult women and men.

Regenerative Technologies Corporation feels a responsibility to provide as much information as possible to potential users. The following listing of precautions was developed from expertise with the technology and potential physiological effects.

The use of the Juvent 1000N DMT Platform has not been evaluated on the following conditions. Individuals with the following conditions or implants should consult their physician before using the medical device:

- Pregnancy (***never use if you know or suspect you are pregnant***)
- Congestive heart failure
- Past history of deep vein thrombosis and/or pulmonary embolism
- History of thrombophlebitis within 5 years
- Sensitivity to motion sickness
- Known retinal conditions (eye)
- Joint implants
- Pacemakers and implantable cardioverter defibrillators (ICDs)
- Treatment/surgery for spinal conditions

NOTE: Some people who are used to very little physical exercise or activity may initially feel tired after the recommended 20-minute treatment. If this occurs after your first treatment, you may wish to start with a shorter treatment time (for example, two minutes or five minutes) and/or rest a day in-between treatments. You can gradually increase the time you use the device over several weeks, until you can comfortably use it for the 20-minute treatment times.